Read PD 1

```
#include "Main.h"
void Read PD (void)
{
   SetDigitalOutput (11, 1); // close shutter, clear film, and increment counter
   SetDigitalOutput (12, 1); // initialize conter value to '0'
   SetDigitalOutput (12,0); // alloww counter to count
   Wait (5); // 5mS wait for things settling down
   SetDigitalOutput (11,0); // open shutter and expose film
   Wait ( expose_time ); // exposer time = 3mS to 8mS
   PD0 = GetAnalogInput (1); // read intensity
   PD1 = expose_and_read();
   PD2 = expose and read();
   PD3 = expose and read();
   PD4 = expose and read();
   PD5 = expose and read();
   PD6 = expose and read();
   PD7 = expose and read();
   PD sum=PD0+PD1+PD2+PD3+PD4+PD5+PD6+PD7;
}
```